

In the back end:
LCD_menu class to
generate the menus

Pressure_Gauge class

Setup
Only executed once

Initialise LCD with
good default values


Calibrate pressure
sensors (take offset)
Calibrate oxygen

Calibrate_motor():
Calibrating the
stepper motor

state_machine()
Inspiration/Expiration

High Priority

Cases 1,2,3 and 4



Medium priority

Cases 1,3

Cases 2,4

Shared variables:
Current State
alarm
pressure
Derived Values

`get_lcd():`
Get input values from
the LCD

`venturi_TidalVolume():`
Integrating the flow rate to
get the tidal volume
measurement per
inspiration

`sound():`
Trigger sound alarm
if necessary

Low Priority

Case 1

Case 2

Case 3

Case 4

`refresh_lcd():`
Refresh the
displayed values to
the lcd

`calc_display_pressure():`
Calculating display
pressure

Change state after time
fixed by doctor / Patient
applies negative pressure
on circuit